




Portable BASIC Balance Master®



Physical Dimensions:

	in (W x D x H)	cm
Dual forceplate	18 x 18 x 2	46 x 46 x 5
Dual forceplate weight	20 lbs.	9 kg.

Electrical Characteristics:

- 100-240 volt / 50-60 Hz / 500 Watt
- Conforms to UL STD 2601-1
- Certified to CAN/CSA STD C22.2 No. 601.1
- Compliant to CE Standards



Options:

- *inVision*™
- *NeuroGames*™
- *D.A.T.a*™

Components:

- BASIC Balance Master® system software
- Dual forceplate
- IBM compatible laptop computer
- Color printer
- Wireless mouse
- Electrical isolation transformer

Accessories Included:

- Rocker board
- Step-up blocks: 4 in (10 cm) and 6 in (15 cm)
- Leveling block: 2 in (5 cm)
- Heel/toe wedges: 6° and 12° A/P
- Inversion/eversion wedges: 3° and 6° M/L
- Foam pad: 18 x 18 x 5 in (46 x 46 x 13 cm)

Specifications subject to change without notice.



Portable BASIC Balance Master®

Standardized Assessment Protocols:

Limits of Stability (LOS) The LOS is an assessment of the voluntary motor system that quantifies impairments in ability to intentionally displace the COG to the patient's stability limits without losing balance. The patient performs the task while viewing a real-time display of their COG position in relation to targets placed at the center of the base of support and at the stability limits. For each of eight directions, the test measures movement reaction time, movement velocity, movement distance, and movement directional control.

modified Clinical Test of Sensory Interaction on Balance (mCTSIB) The mCTSIB is a modification of the original CTSIB that provides information about the patient's ability to maintain postural stability under eyes open firm surface, eyes closed firm surface, eyes open on foam, and eyes closed on foam surface conditions. The mCTSIB enhances an observational test also known as the "Foam and Dome" test by providing an objective measure of patient sway velocity for each of the four task conditions.

Rhythmic Weight Shift (RWS) The RWS quantifies the patient's ability to perform rhythmic movements of their COG from left to right and forward to backward at three distinct paces. During performance of each task, the patient views a real time display of their COG position relative to a target moving at the desired pace and amplitude. For each direction and pace, the RWS measures movement velocity and directional control.

Weight Bearing Squat (WBS) The WBS quantifies the patient's ability to perform squats with the knees flexed at 0°, 30°, 60°, and 90°, while maintaining equal weight on the two legs.

Unilateral Stance (US) The US is a performance test quantifying the patient's ability to maintain postural stability while standing on one leg at a time with the eyes open and closed. The US enhances the observational testing of single leg stance performance by providing an objective measure of patient sway velocity for each of the four task conditions.