

EasyOne - Features and Benefits

The following table lists features and benefits of the EasyOne diagnostic and EasyOne frontline spirometers.

Features	Benefits	Diagnostic	Front-line
Ultrasonic Flow Measurement	Assures fast, accurate, reliable results. No moving parts. No codes to enter	✓	✓
Single patient use spirette	Hygienic mouthpiece reduces risk of cross infection	✓	✓
Alpha Numeric Keypad/Graphic Display	Makes navigation and data entry fast and easy	✓	✓
Graphic Curve Display	Review curves without waiting for printout	✓	✗
Meets ATS 1994 diagnostic standards	Assures accuracy equal to much more expensive systems	✓	✓
Meets NIOSH/OSHA and SSD guidelines	Performs occupational and disability testing to national guidelines	✓	✗
Customizable Configuration	Define test settings and data presentation to suit your needs	✓	✗
Automated Coaching Messages and QC	Helps achieve good patient effort and reliable results	✓	✓
Automated Interpretation	Makes testing faster by eliminating manual computations	✓	✓
Two AA Batteries	NO power cord for true portability. Batteries are easily replaced and readily available	✓	✓
Stores up to 700 test sessions	Review or print results at any time. Facilitates portable testing	✓	✓
Optional 8.5x11 or A4 printout	Single page color printout is easy to read, interpret and store	✓	✓
PC Data Transfer	Easily exports to most popular databases for long term storage and analysis	✓	✓
Toll free Technical Support	Fast accurate answers to all your spirometry questions	✓	✓
3 Year Exchange Warranty	If your spirometer fails, we will replace it overnight, NO charge	✓	✓