


EasyWare v2.8

PC Software for Use with EasyOne Spirometer

EasyWare V2.8 PC Software enables you to use the EasyOne Spirometer in conjunction with a PC running Microsoft Windows 98SecondEdition/ME/2000/XP operating system. The following picture shows a possible test result window:



EasyWare V2.8 features (for details kindly refer to manual):

- Entering, viewing and editing patient data from the PC
- Viewing test results (see screenshot)
- Reviewing of tests including memo field, changing trial acceptability and trial ranking
- Deleting records
- Merging of databases
- Real-time children incentive screen (with Screen Connector only)
- Finding records, i.e. searching for patient name or ID
- Printing and print preview of protocols
- Direct printing from EasyOne via PC
- Exporting data to text files
- Storing data from several EasyOne Spirometers
- Reporting EasyOne Spirometer status, ie serial number, firmware version etc
- Changing EasyOne Spirometer configuration
- Setting of EasyWare preferences

EasyWare V2.8 package (Order No. 2040-2) comprises 1 CD-ROM and 1 User Manual - USB cable is always delivered together with the EasyOne Spirometer

EasyWare V2.8 PC Software is available in English, French, German, Italian and Spanish.