

Printer Compatibility List

In general, a large number of printers are supported by EasyOne™. Due to the large number of printers available on the market, ndd cannot test every printer for compatibility. EasyOne™ supports three different principal driver types: Hewlett-Packard (HP) PCL/(PCL with GUI works, but without guarantee), Epson ESC/P2, and Canon Extended Mode (CEM) commands.

Many laser printer manufactures like Brother, Canon, Kyocera etc. have PCL compatible models and work with the HP b&w setting.

Matrix Printers: Many matrix printer manufactures like Epson, Panasonic etc. have ESC/P2 compatible models and work with the EPSON b&w setting. Please check the specifications of the printer for the ESC/P2 mode.

List of tested printers that work using direct connection to EasyOne™ cradle:

Manufacturer	Printer Model (current models highlighted)	Printer Type	EasyOne Setting	Parallel	USB	Note
Apollo	1000	Inkjet	HP b&w	x		
Brother	B4300, MFC8820D, HL5150D, HL5140 HL2070N	Laser	HP b&w	x	x	
Canon	BJC-620, 2000, 4100, 7000, 7100	Inkjet	Canon b&w / color	x		
Canon	BJC-85, 2100, 2110	Inkjet	Canon b&w / color	x	x	
Canon	S400, S450	Inkjet	Canon b&w / color	x	x	1
Canon (Jap)	BJ-M40	Inkjet	Canon b&w / color	x	x	5
Canon	iP2000, iP2200 is NOT working!, iP4200, iP4300, iP5200, iP90 (portable), iP6100	Inkjet	Canon b&w / color		x	7
Epson	Stylus Color 2, 600, 680, 740, 900	Inkjet	Epson b&w / color	x		
Epson	C60, C62, C64 it NOT working!	Inkjet	Epson b&w / color	x	x	
Epson	C70, C80	Inkjet	Epson b&w	x	x	3
Epson	Stylus Photo	Inkjet	Epson b&w / color	x		
Epson	Stylus Photo 750	Inkjet	Epson b&w / color	x	x	
Epson	LQ 570	Matrix	Epson b&w	x		2
Epson (Jap)	PM780C, PM800C, PM880C, PM950C, PM2200C, PM3000C, PM-G800, PXV700, PX-G900	Inkjet	Epson b&w / color		x	5
HP	LaserJet 4L, 4M, 4MP, 6L, 4050, 4250	Laser	HP b&w	x		
HP	LaserJet 1220, 1200 series, 1150, 1300, 2300 series	Laser	HP b&w		x	4
HP	OfficeJet G85	Inkjet	HP b&w / color	x	x	6
HP	DeskJet 350C, Business inkjet 2230, 2280, 2600	Inkjet	HP b&w / color	x		3
HP	DeskJet 340, 350, 500, 510, 520, 540	Inkjet	HP b&w / color	x		3
HP	DeskJet 550C, 610C, 640C, 840C, 850C, 895C, 940C, 948C, 990, 1120C	Inkjet	HP b&w / color	x		
HP	DeskJet 450, 930, 3820, 5160, 5550, 5551, 5650, 5652, 6122, 7870, 9300	Inkjet	HP b&w / color	x	x	
HP	460CB (portable) (Firmware >= V2.16!) 812C, 842C, 5150, 5151, 5440 , 5740, 5940 , 6540, D5069	Inkjet	HP b&w / color		x	
Oki	8Pplus, B4350	Laser	HP b&w	x	x	
Samsung	ML-7000	Laser	HP b&w	x		
Nashuatec	C7417n (with HP b&w setting)	Laser	HP b&w		x	3

Notes:

1. If EasyOne firmware version is 1.15 or lower, the page size is reduced to ¼ page.
2. Matrix printers in general are very slow.
3. Works only in b&w mode
4. Works only with USB cradle firmware > V01.00.00.00
5. Japanese models
6. Must be set to high resolution, otherwise curve is missing
7. If EasyOne firmware version is 2.14 or lower, the page size is reduced to 1/8 page

Additional Notes:

- All printers are supported if the printing is done 'via PC' using the EasyWare PC software.
- Many low cost laser printers use the GDI language which is not supported by EasyOne™.
- Some matrix printers are supported but the printout is extremely slow.