

## MA 25 / MA 25e

### Air Conduction Audiometer

The MA 25 is engineered to be a portable, tough, and reliable air conduction audiometer with full range of frequencies.

#### Features At a Glance

- Portable air conduction screener (125-8000 Hz)
- AC adapter or battery operated (3 AA batteries)
- Easy fingertip operation with frequency and hearing level controls
- Continuous, pulse, warble tones
- Exceptionally lightweight, only 3.5 lbs
- Soft carrying case included

#### Upgrade Option: MA 25e

- Automated Hughson-Westlake
- Results recalled from internal memory following test
- PC compatible with SNAP software
- Includes patient response switch & USB cable

#### Accurate, Dependable Engineering

Take the MA 25 anywhere with confidence. User-friendly knobs make adjusting frequency and hearing levels easy. Simple push-button operation for right/left, pulse tone and warble.

#### Sleek, Small Design

The super-lightweight MA 25 weighs in at only 3.5 lbs. Small enough to hold in the palm of your hand! Tucked into its carrying case (included), the MA 25 is one of our smallest, lightest portables.



## Technical Data

**Standards:**

Audiometer	EN 60645-1/ANSI S3.6, Type 4
Safety	EN 60601-1
EMC	EN 60601-1-2

**Inputs/Outputs** Patient response switch USB Power Supply

**Rise/Fall Times** ~35 msec

**Distortion** 0.5% typical, 2.5% maximum

**Crosstalk** -100 dB minimum

**On/Off Ratio** Greater than 80 dB

**Frequency Modulation** ±5% frequency modulating rate at 5Hz

**Frequency Range** 125 to 8000 Hz

**Hearing Level Range** -10 to 100 dB HL in 1 or 5 dB steps

**Power Supply** 3 AA batteries/External Power Supply  
5V DC 100, 100 – 240V AC

**Dimensions** 8.9" W x 7.1" D x 2.2" H

**Weight** 3.5 lbs  
(with accessories and carrying case)

**Stimulus** Continuous, Pulse, Warble

**STANDARD ACCESSORIES**

- DD45 audiometric headset with HB7 headband
- Audiogram pad
- Operation manual
- 3 AA batteries
- Carrying case
- AC power adapter

**OPTIONAL ACCESSORIES**

- Patient response switch
- Audiocup headset

**CPT CODES**

92551: Screening, pure tone

92552: Threshold, air

Specifications subject to change without notice.

**Accessories**

Standard Accessories



Carrying Case



Audiogram Pad



AC Adapter



**MAICO Diagnostics**  
10393 West 70th St, Eden Prairie, MN 55344  
Tel: (888) 941-4201 | Fax: (952) 903-4100  
e-mail: info@maico-diagnostics.com | Web: www.maico-diagnostics.com

**ONE YEAR LIMITED WARRANTY:** This warranty is extended to the original purchaser of the instrument, by MAICO, through the Distributor from whom it was purchased. The warranty covers defects in material and workmanship for a period of one year from date of invoice of the instrument to the original purchaser. Accessories which are purchased from MAICO at the same time as the instrument are also under warranty for one year from the date of invoice. For additional information contact your MAICO Distributor.